


---

# Vernier Resistor Board

(Order Code VES-RB)




The Vernier Resistor Board has seven 1 W 5% resistors mounted on a circuit board. Resistor values include 10  $\Omega$ , 15  $\Omega$ , 20  $\Omega$ , 30  $\Omega$ , 39  $\Omega$ , 51  $\Omega$ , and 100  $\Omega$ . Test clip holes on either side of the resistors allow for easy connection to clips on wire leads.

CAUTION: Do not use the Resistor Board on a conducting tabletop.

## What is included in this package?

- Resistor Board
- two sets of wire leads with clips

**NOTE:** Vernier products are designed for educational use. Our products are not designed nor are they recommended for any industrial, medical, or commercial process such as life support, patient diagnosis, control of a manufacturing process, or industrial testing of any kind.



**Measure. Analyze. Learn.™**  
**Vernier Software & Technology**

13979 S.W. Millikan Way • Beaverton, OR 97005-2886  
Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440  
info@vernier.com • www.vernier.com